

Consistency Determination Apparatus



INTRODUCTION

The Consistency Determination Apparatus offers a simple and rapid method of determining the consistency of a pulp suspension in freeness tests. It is more rapid than a Bücher funnel and is suitable for use with most types of pulp.

APPLICATIONS

Pulp

INCLUDED ACCESSORIES

Spreader tool

FEATURES

- Rapid operation
- Easy to use
- Suitable for pulp sample preparation for an automated Laboratory Beater.
- Bench mounted

INSTRUMENT SIZE

Depth: 230 mm (9 in.)
Height: 575 mm (23 in.)
Width: 270 mm (11 in.)

WEIGHT

15 kg (40 lb)

CONNECTIONS

Drainage

Vacuum system to assist drainage (not included)

Electrical: none required

ORDERING INFORMATION

Catalog number 33-35

