


Laboratory Sheet Press



These instruments are used to compress and de-watering lembarcepat made the decision sheet machine. Dioperasikan secara type manual hydraulic plunger pump and pressure types are used to improve gradually

This instrument is used for pressing and de-watering the sheet quickly made on sheet making machine. Hydraulic Type Manually Operated and plunger type pump is employed for raising the pressure gradually.

Feature

Working pressure 3.5 kg/cm²,

Standard Model sheet size dia 165 mm.

Suitable for other size: 160 / 165 / 200 mm diameter.



Laboratory Sheet Press



These instruments are used to compress and de-watering lembarcepat made the decision sheet machine. Dioperasikan secara type manual hydraulic plunger pump and pressure types are used to improve gradually

This instrument is used for pressing and de-watering the sheet quickly made on sheet making machine. Hydraulic Type Manually Operated and plunger type pump is employed for raising the pressure gradually.

Feature

Working pressure 3.5 kg/cm²,

Standard Model sheet size dia 165 mm.

Suitable for other size: 160 / 165 / 200 mm diameter.